# Manual de Uso

# KANMED<sup>®</sup> BABYBED

Art. no: BB1-077/3  
2013-05-17

## **Precaución**

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## 11 EMC guidelines

| Guidance and manufacturer's declaration – electromagnetic emissions   |            |   |
|---|------------|---|
| The BABYBED is intended for use in the electromagnetic environment specified below. The customer or the user of the BABYBED should assure that it is used in such an environment. |            |   |
| Emissions test  | Compliance | Electromagnetic environment - guidance  |
| RF emissions<br>CISPR 11  | Group 1    | The BABYBED uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.  |
| RF emissions<br>CISPR 11  | Class B    | The BABYBED is suitable for use in all establishments including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes. |
| Harmonic emissions<br>IEC 61000-3-2   | Class A    |   |
| Voltage fluctuations/ flicker emissions<br>IEC 61000-3-3  | Complies   |   |

| Guidance and manufacturer's declaration – electromagnetic immunity  |                                  |                                  |   |
|---|----------------------------------|----------------------------------|---|
| The BABYBED is intended for use in the electromagnetic environment specified below. The customer or the user of the BABYBED should assure that it is used in such an environment. |                                  |                                  |   |
| Immunity test   | IEC 60601 test level             | Compliance level                 | Electromagnetic environment - guidance  |
| Electrostatic discharge (ESD)<br>IEC 61000-4-2  | +/- 6 kV contact<br>+/- 8 kV air | +/- 6 kV contact<br>+/- 8 kV air | Floors should be wood, concrete or ceramic tile.<br>If floors are covered with synthetic material, the relative humidity should be at least 30 %. |
| Electrical fast transient / Burst   | +/- 2 kV for power supply lines  | +/- 2 kV for power supply lines  | Mains power quality should be that of a typical commercial or hospital  |

|  |   |   |  |
|--|---|---|--|
| IEC 61000-4-4  | +/- 1 kV for input/output lines   | n/a. for input/output lines   | environment.   |
| Surge<br>IEC 61000-4-5   | +/- 1 kV differential mode<br>+/- 2 kV common mode  | +/- 1 kV differential mode<br>n/a. for common mode  | Mains power quality should be that of a typical commercial or hospital environment.  |
| Voltage dips, short interruptions and voltage variations on power supply input lines<br>IEC 61000-4-11 | <5 % $U_T$<br>(>95 % dip in $U_T$ )<br>for 0,5 cycle<br><br>40 % $U_T$<br>(60 % dip in $U_T$ )<br>for 5 cycles<br><br>70 % $U_T$<br>(30 % dip in $U_T$ )<br>for 25 cycles<br><br><5 % $U_T$<br>(>95 % dip in $U_T$ )<br>for 5 sec | <5 % $U_T$<br>(>95 % dip in $U_T$ )<br>for 0,5 cycle<br><br>40 % $U_T$<br>(60 % dip in $U_T$ )<br>for 5 cycles<br><br>70 % $U_T$<br>(30 % dip in $U_T$ )<br>for 25 cycles<br><br><5 % $U_T$<br>(>95 % dip in $U_T$ )<br>for 5 sec | Mains power quality should be that of a typical commercial or hospital environment. If the user of the [Equipment or System] requires continued operation during power mains interruptions, it is recommended that the [Equipment or System] be powered from an uninterruptible power supply or battery. |
| Power frequency (50/60 Hz) magnetic field<br>IEC 61000-4-8   | 3 A/m   | 3 A/m   | Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment   |

NOTE  $U_T$  is the a.c. mains voltage prior to application of the test level.

Guidance and manufacturer's declaration – electromagnetic immunity

The BABYBED is intended for use in the electromagnetic environment specified below. The customer or the user of the BABYBED should assure that it is used in such an environment.

| Immunity test                 | IEC 60601 test level        | Compliance level | Electromagnetic environment - guidance  |
|-------------------------------|-----------------------------|------------------|---|
|                               |                             |                  | Portable and mobile RF communications equipment should be used no closer to any part of the BABYBED, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.  |
| Conducted RF<br>IEC 61000-4-6 | 3 Vrms<br>150 kHz to 80 MHz | 3 Vrms           | <b>Recommended separation distance</b><br>$d = 1,2\sqrt{P}$   |
| Radiated RF<br>IEC 61000-4-3  | 3 V/m<br>80MHz to 2,5GHz    | 3 V/m            | $d = 1,2\sqrt{P}$ 80 MHz to 800 MHz<br>$d = 2,3\sqrt{P}$ 800 MHz to 2,5 GHz<br><br>where $P$ is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and $d$ is the recommended separation distance in meters (m).<br><br>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, <sup>a</sup> should be less than the compliance level in each frequency range. <sup>b</sup><br><br>Interference may occur in the vicinity of equipment marked with the following symbol.<br><br> |

NOTE 1 At 80MHz and 800MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflected from structures, objects and people.

<sup>a</sup> Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the BABYBED is used exceeds the applicable RF compliance level above, the BABYBED should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the BABYBED.

<sup>b</sup> Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 10 V/m.

**Recommended separation distances between portable and mobile RF communications equipment and the BABYBED**

The BABYBED is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the BABYBED can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the BABYBED as recommended below, according to the maximum output power of the communications equipment.

| Rated maximum output power of transmitter<br>W | Separation distance according to frequency of transmitter<br>m |  |   |
|--|--|--|---|
|  | 150 kHz to 80 MHz<br>$d = 1,2\sqrt{P}$                         | 80 MHz to 800 MHz<br>$d = 1,2\sqrt{P}$ | 800 MHz to 2.5 GHz<br>$d = 2,3\sqrt{P}$ |
| <b>0.01</b>                                    | <b>0.12</b>  | <b>0.12</b>                            | <b>0.24</b>                             |
| <b>0.1</b>                                     | <b>0.38</b>  | <b>0.38</b>                            | <b>0.73</b>                             |
| <b>1</b>                                       | <b>1,2</b>   | <b>1.2</b>                             | <b>2.3</b>                              |
| <b>10</b>                                      | <b>3,8</b>   | <b>3.8</b>                             | <b>7.3</b>                              |
| <b>100</b>                                     | <b>12</b>  | <b>12</b>                              | <b>23</b>                               |

For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

**Note 1:** At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

**Note 2:** These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

**KANMED<sup>o</sup> BABYBED**

KANMED AB  
Gårdsfogdevägen 18B  
SE-168 66 BROMMA – Stockholm  
Sweden

Telephone +46 8 56 48 06 30  
Telefax +46 8 56 48 06 39

E-Mail: [info@kanmed.se](mailto:info@kanmed.se)  
Home page: [www.kanmed.se](http://www.kanmed.se)

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---