KANMED^o WARMINGCABINET

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7. TECHNICAL SECTION EN PROBLEEMOPLOSSING

7.1 TEMPERATURE CONTROL/CALIBRATION

To be performed yearly by a qualified technician only.

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GENERAL INFORMATION

Power Cable: Check the integrity of the power cable and verify proper earth connection.

The working temperature is regulated by the temperature regulator T1 and the over temperature protection by capillary thermostat T2. When checking the calibration of T1 and T2, their value shall be compared to the value of precision thermometer with an air sensor. The sensor T2 can be adjusted through a hole in the bottom of the heating compartment.

Place the sensor of the precision thermometer in the middle of the cabinet. Allow at least 1 hour for proper warming up.

7.1.1 TEMPERATURE REGULATION

INFORMATION ABOUT TEMPERATURE MEASUREMENT AND CALIBRATION

T1 is a powered microprocessor temperature regulator/indicator that has resolution of 1°C.

Corrections can be done through a series of pushing's on the regulator buttons according to the description below.

PROCEDURE

Place the external control sensor as described under GENERAL INFORMATION above.

Start the warming and wait until the temperature is stable (at least 1 hour).

Compare the external control thermometer with the set temperature and if the deviation is bigger than $\pm 1^{\circ}\text{C}$ adjust as follows.

Regulator type No 700-0821

1. Press \uparrow and \downarrow exactly simultaneously (only one beep may be heard, if you don't succeed- .repeat) and keep buttons pressed for 5 sec until the indicator display shows "PA".
2. Press set one time.
3. Press \downarrow within 15s until -19 is displayed.
4. Press set one time.
5. Press \uparrow and \downarrow exactly simultaneously until "SP"(Select Parameter) is displayed
6. To change a parameter press \uparrow or \downarrow button till parameter "CA1" is received.
7. Press SET once and within 2 sec on \uparrow or \downarrow button to enter the noted temperature difference. Press SET again
For example: The display shows 50°C but the external control temperature shows 47°C . The difference is -3°C . To compensate this difference press \downarrow button 3 times. Store the new calibrated value by pressing \uparrow and \downarrow simultaneously (only one beep may be heard, if you don't succeed- repeat) and keep pressed for 5 sec until the indicator display shows the actual temperature.
8. Make a new measurement of the temperature and control that T1 shown value compares to that of the external thermometer.

7.1.2 OVER TEMPERATURE T2 CALIBRATION AND CHECK

The temperature sensor T2 is a capillary thermostat in series with the heating element.

In the Kanmed Cabinet its hysteric is 4°C max. When delivered T2 is set to 5°C above the maximum value of T1. The activation of T2 is identified by a clear click sound as well as when the red over temperature lamp is lit

The procedure described below is a simplified method to check and if necessary adjust T2

PROCEDURE

1. Follow the procedure above, step 1-5.
2. To change a parameter press \uparrow or \downarrow button till parameter "r2" is received. Normally set to 42°C or 70°C . Press "set" and change it to 47 and press set again.
3. Exit the program by pressing the up and down buttons \uparrow and \downarrow until normal view is present.
4. Press Set again and set the temp to 47°C and let the cabinet temperature stabilise.
5. When the temperature is stable adjust the over temperature by turning T2's adjustment screw with a screwdriver (you find it through the hole on the underneath of the shelf). Turn the adjustment screw fully clockwise (+). Now Turn counter clockwise back to activation (one click sounds), turn again clockwise past the activation point in tiny steps (a new click sounds). You have now set the over temperature to about 47°C .
6. Repeat the procedure 1-3 and set back the maximal selectable temperature to 42°C .
7. Check that you don't activate the over temperature indication when running the cabinet in normal operation (with T1 prepared on 42°C) If this should happen you have to repeat the procedure and set T2:s activation point a bit higher.

7.2 TROUBLE SHOOTING

If the warming is not starting, check as follows:

- Power in the wall socket?
- Power Switch on (If power is OK it will show a green light)
- Power cable connected (UK only - plug fuse OK)?
- Cabinet fuses OK?
- Is the fan rotating? It shall start as soon as the power switch is on. If it is not rotating, making a strange noise or rotates slowly it must be changed.
- Check the fan and that the air intake under the heating compartment is not blocked.
- Red lamp lights up now and then. Check the capillary thermostat and the temperature regulating unit.

NOTE: All checking that involves opening the heating compartment must be done by a qualified technician. Ensure that the cables to the heating compartment do not get stuck when you close the heating compartment after service actions

7.3 CHANGING THE MAXIMUM TEMPERATURE, UP OR DOWN.

Follow the procedure described in 7.1.2 and re-adjust T2 so that it is 4-5°C higher than T1.

NOTE: Think carefully about the risks and consequences of setting the max temp limit **higher** than the factory set temperature and about how to make the staff aware that a higher temperature might have been set by someone not authorised to change the temperature.

7.4 ACOUSTIC HIGH TEMP ALARM AND TEMPERATURE LOCK

Please contact Kanmed if you want to activate an acoustic over temperature alarm or lock the possibility to change temperature at all..

8. TECHNICAL DATA

	Large Cabinets		Small Cabinets		Combination
	Universal	Blanket	Universal	Blanket	Fluid and Blanket
Outside measurements Height / Width / Depth	168cm /66cm/ 64cm Height without wheels or feet		95cm / 66cm / 64cm Height without wheels or feet		Same as Large Cabinet
Inner Volume	≈415 litres	≈415 litres	≈190 litres	≈190 litres	≈160 +≈160 litres
Inside measurements Height / Width / Dept	140 cm/56 cm/53 cm		65 cm/56 cm/53 cm		≈59 cm /56 cm/ 53 cm x 2
Height Feet. Height Wheels	Feet 13 -18 cm. Wheels 16cm. On the large Cabinets the front wheels stick out 10 cm and add 2 x 13 cm to the total width.				
Voltages / Power max / Fuses	230VAC ± 10% / ≈900W (1800W double Cabinet) / T6,10A 250V . 115VAC 60Hz is also available. Please contact Kanmed.				
Average power consumption	100 -200W	100 -250W	100 -200W	150-250W	150-250W
Possible number of Basket / Shelves	10	4	5	2	4 in fluid comp. 2 in blanket comp.
Shelf Width x Dept x Height Weight (GE- 41500)	52x50x2 cm / 4 kg	53x49x1,2 cm /2,2 kg	52x50x2 cm / 4kg	53x49x1,2 cm / 2,2 kg	52x50x2 cm / 4kg
Basket size Weight (GE-41600)	52x50X10cm 6kg	Only shelves	52x50x10cm 6 kg	Only shelves	52x50X10 cm / 6kg
Weight without shelves / baskets	≈130 kg	≈130 kg	≈100 kg	≈100 kg	≈140 kg
Doors with magnetic handle	Doors can be hung left or right. Factory standard is left hung. Please specify when ordering. Right hung will increase the delivery time with about 3 weeks.				
Load per shelf/basket	Maximum 20kg				
Warming Capacity	Warms 10 pieces Kanmed Gel GE-455015 from 22°C to 40°C in less than 3 hours.				
Temperature range	35°C to 50°C	35°C to 80°C	35°C to 50°C	35°C to 80°C	36°C to 50°C 36°C to 80°C
Accuracy	± 2°C	± 4°C	± 2°C	± 4°C	±2°C and ±4°C
Factory set maximum temperature	42°C	70°C	42°C	70°C	42°C and 70°C
Other features, information	Made of high quality stainless steel. Polyurethane foam insulation for minimal heat leakage and optimal sound insulation. The Door has double energy saving safety glass. Large and Small Universal Cabinets are equipped with rails for 10/5 shelves or baskets. Height between the shelves and baskets in the Universal Cabinet is ≈13cm. In the Blanket Warming Cabinets the shelves can be moved up or down in steps of 13 cm.				
Temperature regulation	Programmable electronic temperature regulator with display. Temperature can be blocked at a fixed value. A capillary thermostat ensures independent over temperature prevention. A bimetal thermostat is integrated in the heating element.				
Other Kanmed warming products	For anaesthesia. Kanmed Operatherm with Gel Pads and Kanmed WarmCloud. For neonatal use. Kanmed BabyWarmer and Kanmed BabyBed				
CE - marking	According to 89/366/EEC and 93/68/EEC (EN 60601-1-2, EN55011 (1991))				
Expected Lifetime	Kanmed warrants a safe lifetime of 10 years from first day of use. This is under the condition that the Cabinet has been serviced according to the user and or service manual and that the unit has not been modified or changed in any way or for any reason				

9. ACCESSORIES AND SPARE PARTS

Ordering Information					
Article /Ordering Numbers:	Large Cabinet		Small Cabinet		Combination Cabinet
	Universal	Blanket Warming	Universal	Blanket Warming	Fluid and blanket warming
Cabinet with adjustable feet	GE-2350L	GE-2380L	GE-2350S	GE-2380S	GE-2350D
Shelf	GE-41500	-	GE-41500	-	GE-41500
Basket	GE-41600	-	GE-41600	-	GE-41600
Shelf Blanket Cabinet	-	GE-41580	-	GE-41580	GE-41580
Wheel Kit Large Cabinet	GE-41700	GE-41700	-	-	GE-41700
Wheel Kit Small Cabinet	-	-	GE-41900	GE-41900	-


SPARE PARTS

Part no	Position, picture		Description	Quantity
700-0180	4	S2	Mains power switch	1
700-0848	8	M1	Fan	1
700-0862	8	M1	Fan, Only for GE-2350 D (2pcs per cabinet)	1
700-0849	12	HE	Heating Element	1
700-0457	9	G1	Temperature sensor	1
700-0864	9	G1	Temperature sensor in lower compartment GE-2350 D	1
700-0456	10	TR	Transformer (Not valid in GE-2350 D)	1
700-0821	5	T1	Temperature regulator	1
700-0863	5	T1	Temp. regulator GE-2350D Only (2pcs per cabinet)	1
700-0458	6	T2	Thermostat	1
700-0187	7	L1	Red lamp	1
700-0202			Glass for door to cabinet 175 cm high	1
700-0212			Glass for door to cabinet 90 cm high	1
700-0850			Magnetic Handle for door	1
700-0203			Wall mounting bracket	2
700-0669			Door gasket	2m
700-0206			Foot for cabinet	1
			Fuse holder	1
	4		Fuse T6,3A	
	13	F2	Fuse T50mA	
	4		Fuse T 10A Only for 2350 D	
GE-E4055			Pair of rails for cabinets with serial number 4xxx up to 4070	1


Note: Parts and location may vary depending of type of cabinet.

10. CIRCUIT DIAGRAM for cabinets with Art. no. GE-23xx

GE-2350D COMBINATION CABINET:

11. WARRANTY

Kanmed warrants the purchaser that the Warming Cabinet is free from defects in material and workmanship for a period of 24 month from the date of delivery.

The sole obligation of Kanmed with respect to any such defect is limited to the repair with new or re-manufactured parts or, at the discretion of Kanmed, replacement of the equipment or refunding of the purchase price.

This warranty shall not apply if the product has been modified, adjusted or repaired other than by Kanmed or by organisations authorised by Kanmed or modified, adjusted or repaired not in accordance with written instructions provided by Kanmed or if the equipment has been subject to misuse, negligence or accident.

These warranties are made on the condition that prompt notification of a defect is given to Kanmed or its authorised dealers within the warranty period.

Kanmed shall have the sole right to determine whether a defect exists.

Kanmed shall not in any case be liable for special or consequential damages arising from the breach of warranty, breach of contract, negligence or any other legal theory.

When the product have reached end of life, it should be returned to the distributor for recycling in accordance with the EU 2002/96/EC (WEEE) directive if applicable.

Blanket Warming Cabinet

12. EMC COMPATIBILITY STATEMENT

Manufacturer's declaration – electromagnetic emissions	
The KanMed Warmingcabinet GE-2350/GE-2380 is intended for use in the electromagnetic environment specified below. The customer or the user of the KanMed Warmingcabinet GE-2350/GE-2380 should assure that it is used in such an environment.	
Emissions test	Compliance
RF emissions CISPR 11	Group 1
RF emissions CISPR 11	Class B
Harmonic emissions IEC 61000-3-2	Not applicable
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Not applicable
Electromagnetic Environment	The KanMed Warmingcabinet GE-2350/GE-2380 is suitable for use in all establishments, including domestic establishments and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes.

Manufacturer's declaration – electromagnetic immunity	
The KanMed Warmingcabinet GE-2350/GE-2380 is intended for use in the electromagnetic environment specified below. The customer or the user of the KanMed Warmingcabinet GE-2350/GE-2380 should assure that it is used in such an environment.	
IMMUNITY test	IEC 60601 and compliance test level
Electrostatic discharge (ESD) IEC 61000-4-2	<ul style="list-style-type: none"> • ± 6 kV contact • ± 8 kV air
Electrical fast transient/burst IEC 61000-4-4	<ul style="list-style-type: none"> • ± 2 kV for power supply lines
Surge IEC 61000-4-5	<ul style="list-style-type: none"> • ± 1 kV line to line • ± 2 kV line to earth (not applicable)
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<ul style="list-style-type: none"> • $< 5\%$ UT ($> 95\%$ dip in UT) for 0,5 cycle • 40% UT (60% dip in UT) for 5 cycles • 70% UT (30% dip in UT) for 25 cycles • $< 5\%$ UT ($> 95\%$ dip in UT) for 5 s
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	<ul style="list-style-type: none"> • 3 A/m
Conducted RF IEC 61000-4-6	<ul style="list-style-type: none"> • 3 Vrms, 150 kHz to 80 MHz
Radiated RF IEC 61000-4-3	<ul style="list-style-type: none"> • 3 V/m, 80 MHz to 2,5 GHz

KANMED° WARMINGCABINET

Kanmed AB
Gårdsfogdevägen 18B
S-168 67 BROMMA
SWEDEN

Tel +46 (0)8 564 80 630
Fax +46 (0)8 564 80 639

E-Mail: info@kanmed.se
Home page: www.kanmed.se

DISTRIBUTED BY: