

Sweden/Bromma 2014-02-03

### Additional Information Kanmed Water Conditioner BW-50-029

The intended use of the Kanmed Water Conditioner, Art No. BW-50-029, is to prevent bacteria's and algae forming in the Kanmed BabyWarmer Water Mattresses. The 118 ml Water Conditioner is diluted into 4200-10 000 ml of water meaning that the concentration of the various ingredients inside the Water Mattresses will be about 35 to 85 times less than of the contents in the bottle.

Kanmed AB



Lars Karlsson  
Technical Manager