



THE MODERN HOSPITAL BABYBED THAT PROVIDES EFFICIENT ERGONOMICS

The Kanmed BabyBed provides the best care for premature and newborn babies in medium and low intensive care situations. The cosy BabyBed is developed together with neonatal nurses to guarantee a high level of ergonomics. The height of the BabyBed can easily be electrically adjusted and the see-through side walls are designed to fold for easy access to the baby. Soft silicone edges allows you to rest your

arms during treatments. The bed comes with a nicely designed tent to reduce the baby's exposure to light, air condition and sound. Most Kanmed BabyBeds are equipped with the Kanmed BabyWarmer. Together they form a modern care station in Neonatal Intensive Care, Maternity and Delivery wards. Kanmed BabyBed and BabyWarmer are used all over the world.

KANMED BABYBED IS AVAILABLE IN THREE VERSIONS:

(The BabyBeds shown below are equipped with the KANMED BABYWARMER which is ordered separately)



BB-1, KANMED BABYBED

Complete unit includes:

Electrical height adjustment, mattress, tent (model BB-21/3) with tent pole, shelf, power cable, 19 cm high foldable side walls (24 cm height is also available).



**BB-3, KANMED BABYBED
FIXED HEIGHT**

Complete unit includes:

Mattress, tent (model BB-21/3) with tent pole, shelf, 19 cm high foldable side walls (24 cm height is also available).



BB-4, KANMED BABYBED TWINS

Complete unit includes:

Electrical height adjustment, mattress, tent (model BB-22/3) with tent pole, shelf, power cable, 24 cm high foldable side walls. Accommodates 2 complete Kanmed BabyWarmers or 1 Baby Warmer with the Twin Water Mattress and Twin Nest.

ART. NO.	BB-1	BB-3	BB-4
DESCRIPTION AND COLOUR:	Standard bed. Unisex light blue.	Fixed height bed. Unisex light blue.	Twin bed. Unisex light blue.
INCLUDES:	Mattress, tent, tent pole, shelf, 19 cm side walls (24 cm for twin version) power cable (BB-1 & BB-4). Manuals and export packaging.		
OUTER SIZE:	760 mm x 560 mm	760 mm x 560 mm	850 mm x 700 mm
BED SURFACE SIZE/HEIGHT:	665 mm x 445 mm/70 cm-99 cm	665 mm x 445 mm/Fixed at 80 cm	765 mm x 640 mm/70 cm-99 cm
WEIGHT:	About 33 kg	About 25 kg	About 40 kg
WHEEL BASE:	700 mm x 500 mm. Wheels: diameter 125 mm. Front wheels with brakes. Tilt: ≈ 22° C		
MATTRESS:	Polyether with PVC cover. (665 x 445 x 30 mm)	Polyether with PVC cover. (665 x 445 x 30 mm)	Polyether with PVC cover. (765 x 640 x 30 mm)
TENT:	Cotton and polyester. Washable 60° C (all models)		
SHELF:	Min 400 mm x 400 mm. All models. Max load 10 kg		
VOLTAGE:	100-240 Vac, 50/60 Hz	N/A	100-240 Vac, 50/60 Hz
POWER CONSUMPTION:	100 VA	N/A	100 VA
APPLIED STANDARDS: ALL BEDS CARRIES THE CE-MARK	SS EN 60 601. Conforms to MDD.	Conforms to MDD.	SS EN 60 601. Conforms to MDD.


TECHNICAL DATA

KANMED BABYBED ACCESSORIES:



BB-101 DRAWER SYSTEM

Drawer system that is mounted on the back of the bed. Complete with 2 plastic store boxes (13 litres each) and a mounting bracket.



BB-102 BACK SHELF COMPLETE

Vertical poles 25 mm, 2 horizontal rails 10 x 25 mm, 1 top shelf 500 x 300 mm. Including brackets and screws.



BB-TENTS

Different colour tents are available. Please check colour and designs on www.kanmed.se.



700-0743 BOTTLE HOLDER

Stainless steel bracket for holding hand disinfection bottles. Max diameter 70–80 mm.



TENT KIT FOR COTS

Fits into a standard cot. Min. width of cot 34 cm. Can rotate to move tent away from the cot. Is held in place by the mattress.



700-0750 POWER CABLE

Special power cable for powering the Kanmed BabyBed and Kanmed BabyWarmer



700-0829 PULLING HANDLE

To facilitate moving the bed around or for hanging various accessories.

CONCLUSION

The Kanmed BabyWarmer and Kanmed BabyBed allows you to move a premature baby out of the incubator at an early stage and to start the very important mother child bonding using Kangaroo skin to skin care.



KANMED^o BABYWARMER

Kanmed AB

Gårdsfogdevägen 18B | SE-168 66 BROMMA | SWEDEN

Tel +46 8 56 48 06 30 | Fax +46 8 56 48 06 39

info@kanmed.se | www.kanmed.se